

Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days ¹



Control	Lot	
Tri-Level	32406-5K	2024-11-02
Low	32406-51	2024-11-02
Normal	32406-52	2024-11-05
High	32406-53	2024-11-08

MEDONIC

		CONTROL	L
		LOT	32406-51
Medonic M-series/M32		Mean ²	± Range ³
RBC	10 ¹² /L	2.22	0.12
MCV	fL	67.3	7.0
RDW	%	15.8	4.5
Hct	%	14.9	3.0
Plt	10 ⁹ /L	81	15
MPV	fL	11.1	1.5
WBC	10 ⁹ /L	3.4	0.3
Hgb	g/dL	5.5	0.3
MCH	pg	24.8	2.0
MCHC	g/dL	36.9	4.0
Lymph#	10 ⁹ /L	1.5	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.7	0.8
Lymph	%	44.1	7.0
Mid	%	6.0	5.0
Gran	%	49.9	7.0

		CONTROL	N
		LOT	32406-52
		Mean ²	± Range ³
		4.21	0.18
		78.7	7.0
		15.4	4.5
		33.1	3.5
		228	30
		11.7	1.5
		8.7	0.6
		11.6	0.4
		27.6	2.0
		35.0	4.0
		3.9	0.6
		0.6	0.5
		4.2	1.0
		44.5	5.0
		7.3	5.0
		48.2	5.0

		CONTROL	H
		LOT	32406-53
		Mean ²	± Range ³
		5.08	0.22
		93.2	7.0
		13.7	4.5
		47.3	4.0
		518	60
		11.6	1.5
		20.8	1.8
		16.3	0.5
		32.1	2.0
		34.5	4.0
		12.5	1.5
		1.2	1.0
		7.1	2.5
		60.0	5.0
		5.7	5.0
		34.3	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
 - hr Stabilnost otvorene bočice 14 dani (hr)
 - cs Stabilita otevřené ampulky 14 dny (cs)
 - da Åben rør stabilitet 14 dage (da)
 - et Stabiilsus avatud vialis 14 päevade (et)
 - fr Stabilité en flacon ouvert 14 jour (fr)
 - de Stabilität geöffneter Flaschen 14 tage (de)
 - el Σταθιρότητα δειγματος μετά την αποσφράγιση 14 ημέρα (el)
 - it Stabilità della fiala aperta 14 giorni (it)
 - lv Atvērt flakonu stabilitāti 14 dienas (lv)
 - lt Stabilumas atidarius buteliuką 14 dienas (lt)
 - no Åpen rør/glass stabilitet 14 dager (no)
 - pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
 - pt Estabilidade após abertura do frasco 14 dias (pt)
 - ro Termenul de valabilitate al fiolei desfacute 14 zile (ro)
 - ru Стабильность вскрытого флакона 14 дней (ru)
 - sr Proizvod stabilan nakon otvaranja pakovanja 14 dias (sr)
 - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
 - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
 - es Estabilidad de la cápsula abierta 14 días (es)
 - sv Hållbarhet för öppnad flaska 14 dagar (sv)
 - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Sredno
 - Sredstvo
 - Průměr
 - Middelverdi
 - Keskmine
 - Moyenne
 - Mittelwert
 - Μέση
 - Media
 - Nozīmē
 - Vidurkis
 - Middelverdi
 - Średnia
 - Média
 - Medie
 - Среднее
 - Srednja vrednost
 - Priemer
 - Povprečno
 - Media
 - Medelvärde
 - Ortalama

- 3
- Обхват
 - Raspon
 - Rozsah
 - Område
 - Vahemik
 - Intervalle
 - Bereich
 - Ανεπιβρόμενο εύρος α
 - Intervallo
 - Diapazons
 - Intervalas
 - Område
 - Zakres
 - Intervalo
 - Interval
 - диапазон
 - Opseg
 - Rozsah
 - Območje
 - Intervalo
 - Intervall
 - Aralığı



[http://www.medonic.se/
support/downloads](http://www.medonic.se/support/downloads)



25122-4

For Medonic M-series / M32 with barcode reader

M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.:	M-series	Instr.:	M-series	Instr.:	M-series
Lot:	32406-51	Lot:	32406-52	Lot:	32406-53
Prod.:	Con-Diff Low	Prod.:	Con-Diff Norm	Prod.:	Con-Diff High
ID:	3240651+	ID:	3240652+	ID:	3240653+
Exp.:	2-Nov-2024	Exp.:	5-Nov-2024	Exp.:	8-Nov-2024

1:	 Boul e CBAV2	 Boul e CBAV2	 Boul e CBAV2
2:	 3240651+	 3240652+	 3240653+
3:	 24110200001153030000	 24110500002153030000	 24110800003153030000
4:	 10004030800000000000	 10004030100000000000	 10004030800000000000
5:	 22212673700550314930	 42118787701160433135	 50822932701630547340
6:	 08115111152482036940	 22830117152762035040	 51860116153212034540
7:	 03403015040020301708	 08706039060060504210	 20818125150121007125
8:	 15845441700605049970	 15445445500735048250	 13745600500575034350
9:	 11117864830372900898	 11115250915041860444	 11114290979330790113