

Boule Con-5Diff G2

CONTROL

IVD



| Control | Lot |
|-----------|----------|
| Tri-Level | 62107-1K |
| Low | 62107-11 |
| Normal | 62107-12 |
| High | 62107-13 |

Open vial stability:

| |
|------------|
| 2021-12-06 |
| 2021-12-06 |
| 2021-12-06 |
| 2021-12-06 |

14 days¹

Quintus

| WBC | 10 ⁹ /L |
|-------|---------------------|
| NEU# | 10 ⁹ /L |
| LYM# | 10 ⁹ /L |
| MON# | 10 ⁹ /L |
| EO# | 10 ⁹ /L |
| BAS# | 10 ⁹ /L |
| NEU | % |
| LYM | % |
| MON | % |
| EO | % |
| BAS | % |
| RBC | 10 ¹² /L |
| Hgb | g/dl |
| Hct | % |
| MCV | fL |
| MCH | pg |
| MCHC | g/dl |
| RDWcv | % |
| Plt | 10 ⁹ /L |
| MPV | fL |

| LOT L | 62107-11 |
|-------------------|----------------------|
| Mean ² | ± Range ³ |
| 4,1 | 0,5 |
| 2,3 | 0,5 |
| 1,3 | 0,4 |
| 0,3 | 0,3 |
| 0,2 | 0,3 |
| 0,0 | 0,2 |
| 56,4 | 9,0 |
| 32,6 | 9,9 |
| 6,4 | 8,7 |
| 4,4 | 4,5 |
| 0,2 | 1,0 |
| 2,26 | 0,15 |
| 5,2 | 0,5 |
| 15,0 | 2,5 |
| 66,5 | 5,5 |
| 23,0 | 2,8 |
| 34,7 | 3,6 |
| 20,5 | 5,0 |
| 75 | 25 |
| 9,5 | 3,0 |

| LOT N | 62107-12 |
|-------------------|----------------------|
| Mean ² | ± Range ³ |
| 9,0 | 0,8 |
| 4,2 | 1,3 |
| 3,8 | 0,8 |
| 0,6 | 0,7 |
| 0,4 | 0,3 |
| 0,0 | 0,2 |
| 47,3 | 8,0 |
| 41,8 | 9,9 |
| 6,7 | 8,5 |
| 4,0 | 3,0 |
| 0,2 | 1,0 |
| 3,80 | 0,18 |
| 11,9 | 0,6 |
| 32,2 | 3,0 |
| 84,8 | 5,5 |
| 31,3 | 2,8 |
| 37,0 | 3,0 |
| 19,1 | 5,0 |
| 229 | 40 |
| 9,4 | 3,0 |

| LOT H | 62107-13 |
|-------------------|----------------------|
| Mean ² | ± Range ³ |
| 20,6 | 2,0 |
| 7,7 | 3,4 |
| 10,8 | 2,7 |
| 0,9 | 1,9 |
| 1,2 | 0,6 |
| 0,0 | 0,2 |
| 37,3 | 8,0 |
| 52,4 | 9,9 |
| 4,5 | 8,3 |
| 5,6 | 3,6 |
| 0,2 | 1,0 |
| 4,60 | 0,30 |
| 17,6 | 0,7 |
| 45,1 | 3,5 |
| 98,0 | 5,5 |
| 38,3 | 2,8 |
| 39,0 | 3,0 |
| 15,7 | 5,0 |
| 487 | 65 |
| 9,7 | 3,0 |

| | CONTROL | N |
|---|----------|----------|
| i | Target X | Target Y |
| | 99 | 112 |

| | | |
|----|---|--|
| bg | 1 | Трайност на отворен флакон |
| hr | | Stabilnost otvorene bočice |
| cs | | Stabilita otevřené ampulky |
| da | | Åben r r stabilitet |
| et | | Stabiilsus avatud viaalid |
| fr | | Stabilit  en flacon ouvert |
| de | | Stabilit t geoffneter Flaschen |
| el | | Διατ ρηση δειγματος μετ  την αποσφ ραγιση |
| it | | Stabilit  della fiala aperta |
| lv | | Atv rt flakonu stabilit ti |
| lt | | Stabilumas atidarius buteliuka |
| no | | Åpen r r/glass stabilitet |
| pl | | Trwało  po otwarzeniu fiołki |
| pt | | Estabilidade ap s abertura do frasco |
| ro | | Termenul de valabilitate al fiolei desf cut  |
| ru | | Стабильность вскрытого флакона |
| sr | | Proizvod stabilan nakon otvaranja pakovanja |
| sk | | Stabilita otvorenej injek nej liekovky |
| sl | | Stabilnost odprte stekleni ke |
| es | | Estabilidad de la c psula abierta |
| sv | | H llbarhet f r  ppnad flaska |
| tr | | A ık ŐiŐe stabilitesi |

| | | | | |
|----------------|---|------------------|---|---------------------|
| 14 d ni | 2 | Средно | 3 | Обхват |
| 14 dani | | Sredstvo | | Raspon |
| 14 dny | | Pr m r | | Rozsah |
| 14 dage | | Middelv rdi | | Omr de |
| 14 p evade arv | | Keskmine | | Vahemik |
| 14 jour | | Moyenne | | Intervalle |
| 14 tage | | Mittelwert | | Bereich |
| 14 ημ ρα | | Μ ση | | Αναμεν μενο ε ρος α |
| 14 giorni | | Media | | Intervallo |
| 14 dienas | | Noz m  | | Diapazons |
| 14 d. | | Vidurkis | | Intervalas |
| 14 dager | | Middelverdi | | Omr de |
| 14 Liczba dni | | Őrednia | | Zakres |
| 14 dias | | M dia | | Intervalo |
| 14 zile | | Medie | | Interval |
| 14 дней | | Среднее | | диапазон |
| 14 dana | | Srednja vrednost | | Opseg |
| 14 dn  | | Priemer | | Rozsah |
| 14 dni | | Povpre no | | Obmo je |
| 14 dias | | Media | | Intervalo |
| 14 dagar | | Medelv rde | | Intervall |
| 14 g n | | Ortalama | | Araliđı |



www.quintus5part.com

Support/downloads

For Quintus with barcode reader

1. In the Main Menu go to "Calibration QC" / "Enter QC reference".
2. Highlight the first Barcode field, scan in each calibrator/control barcode of the column in sequence (1-9) and click "OK".
3. Repeat steps 1 and 2 for each lot to be entered.

Instr.: Quintus
Lot: 62107-11
Prod.: Con-diff Low
ID: 6210711+
Exp.: 6-Dec-2021

Instr.: Quintus
Lot: 62107-12
Prod.: Con-diff Normal
ID: 6210712+
Exp.: 6-Dec-2021

Instr.: Quintus
Lot: 62107-13
Prod.: Con-diff High
ID: 6210713+
Exp.: 6-Dec-2021

